



P.O. Number:

Date:

2026 Churchill with Positioning Harness Order Form

The Churchill-Pediatric Positioning
Vehicle Restraint System and Accessories

The Churchill is available in four fabric color options. To select a color, simply fill in the space at the end of the part number with the color option abbreviation from the following list.

Sapphire Blue – S
Rapunzel Purple - RP
Total Camo – TC
Beige Naugahyde - N

Billing Address

Contact Name:	
Organization:	
Street Address:	
City, State:	
ZIP Code:	
Phone:	Fax:
E-Mail:	

Shipping Address (If Different)

Contact Name:	
Organization:	
Street Address:	
City, State:	
ZIP Code:	
Phone:	Fax:
E-Mail:	

Qty	Part Number	Description	Unit Price	Line Total
	2000CH15-____	The Churchill with Positioning Harness and 15 Inch Seat Depth – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). NOTE: Add the escape proof Buckle Guard, Chest Clip Guard, and Vehicle Safety Belt Buckle Guard listed below to stop the proficient escape artist from escaping. Also see Note 1 below on how to switch the shoulder belt into locking mode in most vehicles.		
	2000CH18-____	The Churchill with Positioning Harness and 18 Inch Seat Depth – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). NOTE: Add the escape proof Buckle Guard, Chest Clip Guard, and Vehicle Safety Belt Buckle Guard listed below to stop the proficient escape artist from escaping. Also see Note 1 below on how to switch the shoulder belt into locking mode in most vehicles.		

	NOTE	The following three items are designed to stop the “Escape Artist” from escaping. Also see Note 1 below on how to switch the shoulder belt into locking mode in most vehicles.		
	2000BG	Buckle Guard – Stops the “Escape Artist” from opening crotch buckle		
	2000CCG-____	Chest Clip Guard – Stops the “Escape Artist” from opening Chest Clip Guard and from pushing it down		
	8000VSBBG	Vehicle Safety Belt Buckle Guard – Stops the “Escape Artist” from opening the vehicle seat belt. We also advise switching the Retractor to automatic locking retractor mode *****See Note 1 below.		
	2000EUHR-____	EZ-Up Head Rest – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below. The cap is designed to Velcro to the head rest.		
	See note to the right regarding the five cap color options.	Note: There are five options for cap color. The Khaki option is the coolest option for hot temperatures. Simply fill in the space at the end of the part number below with “K” for Khaki or one of the color abbreviations (to match your seat cover) at the top of this form.		
	1000EUCS-____	EZ-Up Cap – Size Small - (Head Circumference: 17.5 – 19.5 inches)		
	1000EUCM-____	EZ-Up Cap – Size Medium - (Head Circumference: 19.5 – 21.5 inches)		
	1000EUCL-____	EZ-Up Cap – Size Large - (Head Circumference: 21 – 23.25 inches)		
	1000EUCC-____	Custom Cap – Please Provide Following Measurements: Head Circumference: _____; Distance from top of forehead over scalp to hairline at back of head _____		
	2000FA-____	Foam Abductor (also known as a pommel)		
	2000HIW-____	Hip Inflection Wedge (for hip flexion)		
	2000HICP-____	Hip Inflection Comfort Pommel (for hip flexion and abduction)		

	2000IC15-____	Quick-Change Incontinent Cover for 15 Inch Seat Depth Base		
	2000IC18-____	Quick-Change Incontinent Cover for 18 Inch Seat Depth Base		
	2000ADD	Adductor Strap		
	2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – (this part is usually only necessary when the vehicle is manufactured prior to 2003, or the Churchill needs to be installed in the front seat. See Note 2 below.		
	NOTE	One of the following positioning vests can be added for simultaneous use with the positioning harness. This would also allow the possibility of using only the positioning vest with the Churchill; or only the positioning harness with the Churchill. If you just want the Churchill with positioning vest, please use the order form on that product page; which will avoid a vest upcharge.		
	2000CVS-____	Small Positioning Vest - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches.		
	2000CVM-____	Medium Positioning Vest - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches.		
	2000CVL-____	Large Positioning Vest - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches.		
	2000CVC-____	Custom Vest – Please provide the following measurements of seated occupant: Chest Circumference: ____; Waist Circumference: ____; Distance from top of thigh to top of shoulder: ____		
	2000VB	Add 4 side release buckles to one of the vests above to make it more difficult to escape the vest.		
			Sub Total**	
Indiana Residents Only; Not Applicable to other States			Sales Tax (7%)	
			TOTAL	

*****The “Churchill with Positioning Harness” includes the following features:**

- Base (15 inch or 18 inch Seat Depth)
- Lower Anchor Attachment Straps
- Support Harness for 5-Point Positioning Harness
- 5- Point Positioning Harness
- Foam Padding
- Fabric Cover
- Tether Extension Assembly
- Belt Positioning Guides

¹ To help contain the “escape artist”, it is important to “switch” the shoulder belt to the locking mode when possible. Many vehicles have a “switchable retractor”. To determine if you have one, slowly pull all of the shoulder belt webbing out of the retractor. Once you pull it all out, begin “feeding” it back in. If you hear a clicking noise, or it remains locked when you gently pull on it, then you have a switchable retractor. Follow the steps in the Churchill’s instructions to “switch” or “lock” the retractor after the occupant is installed in the seat.

² The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH, or the desired seating position doesn’t have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket Kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket Kit will not be necessary. Also, many 2000 – 2002 models also have LATCH (See your vehicle owner’s manual to verify if LATCH is present for these model years). For vehicles without LATCH, the Attachment Bracket Kit can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket Kit cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle’s floor board.

Please call Customer Service at 317-409-0148 with any questions regarding the installation of the Churchill.

Please place orders by one of the following methods:

E-Mail: order@merrittcarseat.com

Phone: 317-409-0148

Fax: 317-893-2567

Thanks for Your Business!

Merritt Manufacturing, Inc.
www.merrittcarseat.com

580 Tracy Rd,
New Whiteland, IN 46184